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食品中黴菌毒素檢驗方法－黃麴毒素之檢驗 Draft Method of Test for Mycotoxin in Foods－Test of Aflatoxins

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Draft Method of Test for Mycotoxin in Foods－Test of Aflatoxins

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2.8.18. DETERMINATION OF AFLATOXIN B₁ IN HERBAL DRUGS

CAUTION: aflatoxins are very toxic and carcinogenic. Perform manipulations under an extraction hood whenever possible. Take particular precautions, such as use of a glove box, when toxins are in dry form because of their electrostatic properties and the tendency to disperse through the working areas. Decontamination procedures for laboratory wastes of aflatoxins were developed by the International Agency for Research on Cancer (IARC).

- Aflatoxins are naturally occurring mycotoxins produced mainly by *Aspergillus flavus* and *Aspergillus parasiticus*. These fungi are common and widespread in nature and are most often found when certain grains are grown under conditions of stress such as drought. The mould occurs in soil, decaying vegetation, hay, and grains undergoing microbial spoilage, and invades all types of organic substrates whenever and wherever the conditions are favourable for its growth. Favourable conditions include high moisture content and high temperature. At least 13 different types of aflatoxin are produced in nature and most of these are known to be highly toxic and carcinogenic. Aflatoxin B₁ is considered the most toxic. Herbal drugs that are subject to contamination by aflatoxins are tested by a validated method.

Unless otherwise indicated in the monograph, herbal drugs contain not more than 2 µg/kg of aflatoxin B₁. The competent authority may also require compliance with a limit of 4 µg/kg for the sum of aflatoxins B₁, B₂, G₁ and G₂.

Aflatoxin	Limit
B ₁	2 µg/kg
sum of B ₁ , B ₂ , G ₁ , G ₂	4 µg/kg

The method described below is cited as an example of a method that has been shown to be suitable for devil's claw root, ginger and senna pods. Its suitability for other herbal drugs must be demonstrated or another validated method used.

METHOD

Liquid chromatography (2.2.29).

Aflatoxins are subject to light degradation. Carry out the determination protected from daylight by using UV-absorbing foil on windows in combination with subdued light, or curtains or blinds in combination with artificial light (fluorescent tubes are acceptable). Protect aflatoxin-containing solutions from daylight.

Rinse glassware before use with a 10 per cent V/V solution of *sulfuric acid R* and then rinse carefully with *distilled water R* until no more acid is present.

Test solution. Use an immunoaffinity column containing antibodies against aflatoxin B₁ with a capacity of not less than 100 ng of aflatoxin B₁ and which gives a recovery of not less than 80 per cent when a solution of 5 ng of aflatoxin B₁ in a mixture of 12.5 mL of *methanol R* and 87.5 mL of *water R* is passed through. Allow the immunoaffinity column to reach room temperature. To 5.00 g of the powdered drug (500) (2.9.12) add 100 mL of a mixture of 30 volumes of *water R* and 70 volumes of *methanol R* and extract by sonication for 30 min. Filter through folded filter paper. Pipette 10.0 mL of the clear

Aflatoxin B₁ standard solution
250 µg/mL
at room temperature
each disc

S

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Calibration solution
B₁ standard solution
1.8 µg/mL
for liquid chromatography
must be used

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Comparison of performance of VICAM affinity columns

Brand name	Analyte	Gel bead [mL]	Column capacity [ng]	Spike level for recovery ng	Spike conditions for recovery	Recovery %
AflaTest®	B1, B2, G1, G2, M1, M2	0.250	300	100 total	Spiked sample extract	B1, G1 ≥ 90%; B2 ≥ 85%; G2 ≥ 80%
AflaTest® WB	B1, B2, G1, G2, M1, M2	0.300	400	200 total	Spiked sample extract	B1, G1 ≥ 90%; B2 ≥ 85%; G2 ≥ 80%
AflaTest® WB SR	B1, B2, G1, G2, M1, M2	0.375	500	400 total	Spiked sample extract	B1, G1 ≥ 90%; B2 ≥ 85%; G2 ≥ 80%
Afla M1™ HPLC	M1, M2	1.000	150	2.5	spiked skim milk	M1 ≥ 85%
OchraTest™	OTA	0.200	100	20	spiked extract	B₁, G₁ ≥ 90% B₂ ≥ 85 % G₂ ≥ 80 %
OchraTest™ WB	OTA	0.250	300	20	spiked extract	
ZeralaTest™	ZON	0.350	1500	1500	spiked extract	
ZeralaTest™ WB	ZON	0.350	1500	1500	spiked extract	
DONtest™	DON	0.400	1250	25	spiked water	≥ 85% T-2
DONtest™ WB	DON	0.450	1250	25	spiked water	
T-2test™ HPLC	T2, HT2	0.300	1500 T-2	1500 T-2	spiked PBS	≥ 85% of each
T-2/H-T2™ HPLC	T2, HT2	0.600	2000 total	1400 total	spiked PBS	≥ 85%
FumoniTest™	FB1, FB2 & FB3	0.300	10000	7000	spiked PBS	≥ 85%
FumoniTest™ WB	FB1, FB2 & FB3	0.350	10000	7000	spiked PBS	≥ 85%
AflaOchra HPLC™	B1, B2, G1 & G2	0.500	100	100	spiked sample extract	> 85%
	OTA		100	100		> 85%
AOZ HPLC™	B1, B2, G1, G2	0.700	100	100	spiked sample extract	> 85%
	OTA		100	100		> 85%
	ZON		1000	50		> 85%
CitriTest® HPLC	Citrinin	0.350	20	20	spiked phosphoric acid solution	> 80%

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